

Rapid Detection of *Salmonella* Using Real-Time PCR Assay in Meat, Poultry, and Whole Liquid Egg Enriched with an Improved Culture Broth

Marie-Christine Ety¹, Marie Goreth Nicizanye¹, Smina Messaoudene¹, Anna Yattara¹, Anne Helmer¹, Alex Charbonneau¹, Sergiy Olishevskyy¹

¹FoodChek Laboratories Inc., St-Hyacinthe, Quebec, Canada


Introduction

Although rapid testing of *Salmonella* in foods has advanced over the past two decades, the development and implementation of rapid and user-friendly, plus accurate methods, remains to be an important task for the food industry.


The objective of this study was aimed on improving Actero™ *Salmonella* Enrichment Media and optimizing the enrichment conditions for the rapid detection of *Salmonella* in foods using a RT-PCR method.

In vitro Culture Studies

Uninjured *Salmonella*



Heat-Injured *Salmonella*



Culture Conditions

- Initial Inoculum – 20 CFU/well
- Media Volume – 200 µL
- Temperature – 35°C
- Time – 18 hours


- Actero™ *Salmonella* 1.1
- Actero™ *Salmonella* 1.2
- BPW
- UPB

Table 1. Growth Kinetics Data


Bacterium	Kinetics Parameter	Actero™ <i>Salmonella</i> 1.1	Actero™ <i>Salmonella</i> 1.2	Buffered Peptone Water (BPW)	Universal Pre-Enrichment Broth (UPB)
Uninjured <i>Salmonella</i>	LPD	6.77 ± 0.23	6.13 ± 0.31*	6.62 ± 0.04	6.93 ± 0.19
	GR	0.66 ± 0.03	0.73 ± 0.03	0.85 ± 0.09**	0.71 ± 0.07
Heat-injured <i>Salmonella</i>	LPD	8.02 ± 0.51	7.35 ± 0.29	8.46 ± 0.33	8.28 ± 0.58
	GR	0.67 ± 0.06	0.71 ± 0.07	0.73 ± 0.06	0.64 ± 0.07

Notes: LPD – Lag Phase Duration; GR – Growth Rate.
* P<0.05; ** P<0.01 as compared to the Actero™ *Salmonella* 1.1 values.

Assay Principle



AOAC Validation Design



Matrix Studies

Table 2. Food Matrices and Inoculating Microorganisms

Food Matrix	Strains (serogroup)	Origin	MPN, CFU/sample	Stress	Normal flora, CFU/sample
Raw Ground Beef	S. Carmel (J)	N/D	0.7	N/A	6.0x10 ⁴ – 8.15x10 ⁵
Chicken Carcass Rinse	S. Enteritidis (D1)	Chicken Wings	1.4	N/A	1x10 ⁷
Whole Liquid Egg	S. Heidelberg (B)	Chicken Legs	0.7	Heat	N/A

Notes: N/A – Not Applicable.

Table 3. Food sampling and Enrichment Conditions

Food Matrix	Sample Size	Medium Volume	Malachite green	Homogenization	Incubation temperature	Incubation Time	Dilution for PCR
Raw Ground Beef*	375 g	1125 mL (1/4)	25 mg/L	60 sec	39°C	20 hrs	1:25
Chicken Carcass Rinse	30 mL	30 mL (1/2)	20 mg/L	Manually	35°C	16 hrs	1:25
Whole Liquid Egg*	100 g	300 mL (1/4)	N/A	Manually	35°C	18 hrs	1:25

Notes: * pH was adjusted to 6.8 using 10 N NaOH; N/A – Not Applicable.

Table 3. Detection of *Salmonella* in Meat, Poultry, and Whole Liquid Egg Using BAX® System Assay - Method Verification Study

Food Matrix	Media	N	Candidate Method (C)		Reference Method (R)	
			RT-PCR POD	Culture Method POD	MLG 4.09 POD	dPOD (C,R) (95%CI)
Raw Ground Beef	Actero <i>Salmonella</i> 1.1	30	0.60	0.60	0.45	0.15 (-0.15; 0.41)
	Actero <i>Salmonella</i> 1.2	30	0.75	0.75	0.65	0.10 (-0.18; 0.36)
Chicken Carcass Rinse	Actero <i>Salmonella</i> 1.1	30	0.85	0.85	0.75	0.10 (-0.15; 0.34)
	Actero <i>Salmonella</i> 1.2	30	0.80	0.80	0.75	0.05 (-0.21; 0.30)
Whole Liquid Egg	Actero <i>Salmonella</i> 1.1	30	0.60	0.60	0.55	0.05 (-0.24; 0.33)
	Actero <i>Salmonella</i> 1.2	30	0.40	0.40	0.55	-0.15 (-0.41; 0.15)

Notes: POD - Probability of Detection.

Conclusions

- ✓ The candidate method using Actero™ *Salmonella* Enrichment Media for recovery of *Salmonella* in meat, poultry and eggs followed by the detection using BAX® System RT-PCR Assay showed equivalent performance in comparison with the USDA-FSIS MLG 4.09 reference method.
- ✓ The incorporation of specific selective growth substrates and stimulating co-factors, previously used as liquid supplements into a powdered formulation of Actero™ *Salmonella* Enrichment Media, could improve the capacity of *Salmonella* to grow in the broth.

