

Actero™ EZ-Media Liquid Bag



“Always Easy. Always Right.”

For all your food and environmental
pathogen testing needs.



foodchexsystems.com

Increases Safety

Contains sterile enrichment medium that eliminates autoclaving or hot glassware.

Increases Throughput

Ready-to-use liquid medium that is quick to dispense through the supplied filter.

Reduces Human Error

Contains pre-mixed liquid medium – no media preparation is required.

Reduces Time

Eliminates measuring, heating or standard consumables.

High Quality

Ingredients from Europe and U.S.A.

Product List

Pathogen	FoodChek Cat. #	Product Name	Format	Description
Listeria	FCM-203	ACTERO™ ELITE LISTERIA ENRICHMENT MEDIA – EZ-MEDIA LIQUID BAG	3x3L	The patented Actero™ Elite Listeria Enrichment Media is a selective medium, specifically optimized for single-step recovery and enrichment of <i>Listeria</i> spp. from food and environmental samples. The addition of ferric ammonium citrate allows for the improved recovery of <i>Listeria</i> spp. A box of 3x3L or 2x5L liquid Actero™ Elite Listeria Media.
	FCM-204		2x5L	
	FCM-242	ACTERO™ DEMI FRASER – EZ-MEDIA LIQUID BAG	3x3L	Demi-Fraser Media is a modification of Fraser Broth Base, this media was developed by Fraser and Sperber for the rapid detection of <i>Listeria</i> from food and environmental samples. The addition of ferric ammonium citrate allows for the improved recovery of <i>Listeria</i> spp. A box of 3x3L or 2x5L liquid Actero™ Demi Fraser Media.
	FCM-243		2x5L	
Salmonella	FCM-239	ACTERO™ BUFFERED PEPTONE WATER ISO – EZ-MEDIA LIQUID BAG	3x3L	Buffered Peptone Water ISO is used for pre-enriching damaged <i>Salmonella</i> spp. from various food sources. A box of 3x3L or 2x5L liquid Actero™ ISO BPW Media.
	FCM-241		2x5L	



foodcheksystems.com

“Always Easy. Always Right.”

You can always depend on the Actero™ brand of enrichment media products to deliver high quality and reliable performance with time and cost savings.

FoodChek[™]
Food safety, **simplified.**