

Actero™ Oxford Listeria Agar Product Information

Catalogue No.	Description
FCM-084	Actero™ Oxford Listeria Agar (500 G)
FCM-083	Actero™ Oxford Listeria Agar (2 KG)
FCM-082	Actero™ Oxford Listeria Agar (10 KG)

INTENDED USE:

With the addition of antibiotics, Oxford Listeria Agar is used for the selective isolation of *Listeria monocytogenes* from food. Typically, *Listeria monocytogenes* is visible after 24 hours but continued incubation may allow for the detection of slower-growing strains. *Listeria* hydrolyzes esculin to produce black zones around the colonies while Lithium chloride and antibiotics are incorporated in the medium for selectivity.

Formula* per Liter:

Columbia Blood Agar Base	39.0g
Esculin.....	1.0g
Ferric Ammonium Citrate	0.5g
Lithium Chloride.....	15.0g
Agar.....	2.0g

Final pH: 7.2 ± 0.2 at 25°C

* Grams per liter may be adjusted or formula supplemented to obtain the desired performance.

PREPARATION:

Mix 57.5 grams of the medium in one liter of purified water until evenly dispersed. Heat with repeated stirring and bring to a boil to dissolve completely. Distribute and autoclave for 15 minutes at 121°C.

TERMS AND CONDITIONS:

Salus Scientific Inc. makes no representations and warranties concerning its products other than those stated herein. All Product(s) delivered hereunder by Salus Scientific Inc., its affiliates or any other person on its behalf shall, at the time of delivery, be manufactured to meet Salus Scientific Inc.'s specifications and all applicable laws. All other terms, conditions and warranties, including any warranty of merchantability, quality, fitness or suitability for a particular or intended purpose, implied by common law or statute, (implied warranties) are expressly excluded.

FOR FURTHER INFORMATION, PLEASE CONTACT:

Salus Scientific Inc.
Suite 450, 1414 – 8 St. S.W.
Calgary, Alberta, T2R 1J6
Canada
Tel: 1-877-298-0208

QUALITY CONTROL SPECIFICATIONS:

1. The powder is homogeneous, free flowing, and beige.
2. Visually the prepared is light to medium beige and slightly hazy.
3. Expected cultural response after 24-48 hours at 35°C. Media prepared with Oxford Antimicrobial Supplement and Modified Oxford Antimicrobial Supplement.

Organism	Oxford	Modified Oxford
<i>Escherichia coli</i> ATCC® 25922	Inhibition	Inhibition
<i>Listeria monocytogenes</i> ATCC® 7644	Growth, Black colonies	Growth, Black colonies
<i>Listeria monocytogenes</i> ATCC® 19114	Growth, Black colonies	Growth, Black colonies
<i>Staphylococcus aureus</i> ATCC® 25923	Suppressed	Inhibition

Storage Instructions:

Store the sealed bottle containing the dehydrated medium at 2 to 30°C. Once opened and recapped, place the container in a low humidity environment at the same storage temperature. Protect it from moisture and light. The dehydrated medium should be discarded if it is not free flowing, or if the color has changed from the original beige color.